



# Affinity Pixie™

OXYGENATION SYSTEM  
WITH BALANCE™ BIOSURFACE



Performance  
and  
Possibilities  
for Your Small  
Patients

Innovating for life.


INTRODUCING THE

# Affinity Pixie™

OXYGENATION SYSTEM

Our addition to the Affinity® family

The comprehensive system design of the oxygenator, cardiotomy/venous reservoir and holder provides ease of set-up, ease of use and minimized prime volumes.



## Balance™ Biosurface


**Thromboresistance and enhanced blood compatibility**

- A hydrophilic biosurface option without heparin
- Applied to primary blood contacting surfaces of the oxygenator and cardiotomy/venous reservoir

## The Affinity Pixie™

Oxygenation System delivers performance and versatility for neonates, infants and small children requiring cardiopulmonary bypass at flow rates up to 2.0 L/min.

This addition to the Affinity family is part of Medtronic's commitment to providing more options for pediatric patients undergoing cardiopulmonary bypass.






# Oxygenator

Advanced Affinity® oxygenator technology in a compact, low-prime design

- 2.0 L/min maximum flow rate
- 48 mL prime
- Plasma resistant polypropylene microporous fiber membrane
- Fiber bundle wound directly over PET heat exchanger for reduced prime volume
- Indicated for transfer of volatile anesthetics isoflurane, sevoflurane, desflurane and enflurane
- Device and accessories do not contain DEHP



Blood contact:  
blood side of  
pre-membrane  
purge port

Water contact:  
water pathway  
above heat  
exchanger

## Electrical shunt


This safety feature dissipates electrostatic charge that may develop inside the oxygenator's blood pathway due to rotation of the roller pump heads.

Electrostatic charge is shunted from the blood path to the water path where it is discharged through the heater/cooler unit's electrical grounding system.


# Oxygenator Performance Data

Data on file – ISO Standard Conditions


## O<sub>2</sub> Transfer Rate




## CO<sub>2</sub> Transfer Rate



## Heat Exchange Effectiveness



## Blood Side Pressure Drop



The oxygenator's proprietary graduated fiber bundle density technology together with radial flow path design provide enhanced gas transfer, low blood-side pressure drop, decreased prime volume, uniform blood flow distribution and avoidance of areas of stasis.







# Cardiotomy/Venous Reservoir

**Ergonomic, easy to use design**  
**Minimized blood contact with antifoam**

- 2.0 L/min maximum flow rate
- 1200 mL volume capacity
- 20 mL minimum operating volume without vortex formation
- Independent rotation of venous inlet and cardiotomy turret
- Port locations and spacing allow easy tubing attachment
- Vacuum assisted venous drainage-ready design
- Device and accessories do not contain DEHP

**Minimized blood contact with antifoam**



## FEATURING THE Affinity Orbit™ Holder System

**Quick set-up. Custom device positioning and orientation**

- Quick and easy device attachment
- Adapts to any pump set-up
- Versatile device positioning and port orientation
  - Holder arm, oxygenator arm and cardiotomy/venous reservoir pivot for custom device positioning
  - Oxygenator rotates 360° for desired port orientation
- Facilitates reductions in circuit tubing length and associated prime volume
- Designed for use with future Medtronic perfusion technologies



Custom device positioning and orientation



### Affinity® Manifold Holder

- Holds a standard flat-plate stopcock manifold in either horizontal or vertical orientation
- C-clamp for easy and stable attachment to both IV poles and larger diameter pump masts

Even greater  
**positioning flexibility**  
to adapt to your preferred circuit setup.



# Affinity Pixie™ Oxygenation System

## ORDERING INFORMATION

Oxygenation Systems		
Model Number	Description	Units/case
BBP211	Hollow Fiber Oxygenator with Balance™ Biosurface	4
BBP241	Hollow Fiber Oxygenator and Cardiotomy/Venous Reservoir with Balance™ Biosurface	4
Accessories and Holders		
Model Number	Description	Units/case
ATP210	Affinity® Temperature Probe	1
AUH2093	Affinity Orbit™ Holder System	1
AMH2014	Affinity® Manifold Holder	1

\*Technology licensed under agreement from BioInteractions, Limited, United Kingdom. Balance™ is a trademark of Medtronic, Inc.

For information on other Medtronic technologies for extracorporeal circulation, blood processing and diagnostics, visit: [www.perfusion.medtronic.com](http://www.perfusion.medtronic.com)

**Related products:** Contact your Medtronic Sales Representative for information on additional Medtronic technologies for pediatric cardiopulmonary bypass, including arterial filters, cannulae products and custom perfusion systems incorporating Medtronic pediatric CPB devices

The Affinity Pixie™ Oxygenation System is indicated for use in an extracorporeal circulation circuit during cardiopulmonary bypass procedures up to 6 hours in duration. A strict anticoagulation protocol should be followed and anticoagulation should be routinely monitored during all procedures. The benefits of extracorporeal support must be weighed against the risk of systematic anticoagulation and must be assessed by the prescribing physician.

For a listing of indications, contraindications, precautions and warnings, please refer to the Instructions for Use which accompany each product.

## SPECIFICATIONS

Oxygenator	
Membrane type	Microporous polypropylene hollow fiber
Membrane surface area	0.67 m <sup>2</sup>
Heat exchanger	Polyethylene Terephthalate (PET)
Static priming volume	48 mL
Recommended blood flow range	0.1 – 2.0 L/min
Maximum water pressure	1550 mmHg (206 kPa)
Maximum blood pressure	750 mmHg (100 kPa)
Venous blood inlet	1/4 in (0.6 cm)
Arterial blood outlet	1/4 in (0.6 cm)
Arterial sample port	Female luer
Recirculation port	Female luer
Pre-membrane purge port with one-way valve	1/16 in (0.16 cm), male luer
Biosurface	Balance™ Biosurface
Accessories	
Tubing adapter	3/16 in (0.5 cm) (2)
Recirculation line	3/16 in (0.5 cm) (1)
Cardiotomy/Venous Reservoir	
Reservoir volume capacity	1200 mL
Recommended blood flow range	0.1 – 2.0 L/min
Maximum cardiotomy flow rate	2.0 L/min
Minimum operating level	20 mL
Cardiotomy filter	30 µm nominal polyester depth filter
Venous inlet screen	64 µm
Venous return inlet, rotatable	1/4 in (0.6 cm)
Venous reservoir outlet	1/4 in (0.6 cm)
Venous inlet luer locks	2 luer locks
Cardiotomy inlet, step-up	3/16 in (0.5 cm) to 1/4 in (0.6 cm) (3)
Cardiotomy inlet	1/4 in (0.6 cm) (1)
Filtered luer locks to cardiotomy filter	4 luer locks
Nonfiltered luer locks	2 luer locks
Filtered quick prime port	1/4 in (0.6 cm)
Nonbarbed vent port	1/4 in (0.6 cm)
Maximum rated pressure	+ 20 mmHg/ - 100 mmHg
Positive/negative pressure relief valve	< 5 mmHg positive/ 60 mmHg vacuum
Biosurface	Balance™ Biosurface
Accessories	
Tubing adapter	3/16 in (0.5 cm) (2)
Tubing adapter	3/8 in (1.0 cm) (2)
Flexible luer lock adapter	2 adapters
Sampling manifold	1 manifold assembly

### World Headquarters

Medtronic, Inc.  
710 Medtronic Parkway  
Minneapolis, MN 55432-5604  
USA  
Tel: (763) 514-4000  
Fax: (763) 514-4879

### Europe

Medtronic International Trading Sàrl  
Route du Molliat 31  
Case postale  
CH-1131 Tolochenaz  
Switzerland  
Tel: 41.21.802.7000  
Fax: 41.21.802.7900

### Canada

Medtronic of Canada Ltd.  
99 Hereford Street  
Brampton, Ontario L6Y 0R3  
Canada  
Tel: (905) 460-3800  
Fax: (905) 826-6620  
Toll-free: 1 (800) 268-5346

### Asia Pacific

Medtronic International, Ltd.  
16/F Manulife Plaza  
The Lee Gardens, 33 Hysan Avenue  
Causeway Bay  
Hong Kong  
Tel: (852) 2891 4456  
Fax: (852) 2891 6830  
[enquiryap@medtronic.com](mailto:enquiryap@medtronic.com)

### Latin America

Medtronic USA, Inc.  
Doral Corporate Center II  
3750 NW 87th Avenue Suite 700  
Miami, FL 33178  
USA  
Tel: (305) 500-9328  
Fax: (786) 709-4244



**Medtronic**