

Amphirion™ Deep

Infrapopliteal PTA Balloon Catheter OTW 0.014"

Easy access to the extremities*

- Tip profile (0.017")
- LFC hydrophilic coating
- Three-step shaft design to optimize tracking through distal vessels
- Reinforced proximal shaft design for strong pushability
- 4 F compatible in all sizes

Balloon sizes to accommodate your needs*

- Size range 1.5–4.0 mm in diameter and 20–210 mm in length
- Long balloons (up to 210 mm) are suited for treatment of extremely diffuse lesions

Conformable balloon material*

- Proprietary polymer blend provides wonderful conformability

Tapered balloon*


- The 210 mm balloon tapers by 0.5 mm to respect the arterial anatomy

TECHNICAL SPECIFICATIONS

Catheter design	Over the Wire (OTW)
Balloon marker	1/2 swaged, Platinum Iridium
Shaft diameter	Proximal 3.9 F Middle 3.3 F Distal 2.8 F
Introducer sheath compatibility	4 F
Guidewire compatibility	0.014"
Nominal pressure	7 atm



Image courtesy of
Dr. Marco Manzi, Italy.



Below the Ankle



* Bench Test Data on file with Medtronic, Inc.
Test data not indicative of clinical performance

Amphirion™ Deep

Infrapopliteal PTA Balloon Catheter OTW 0.014"

ORDER INFORMATION

Product Catalogue Number		Balloon Diameter (mm)	Balloon Length (mm)	Recommended Introducer Sheath (F)	Nominal Pressure (BAR)	RBP (BAR)
OTW Usable Length 120 cm	OTW Usable Length 150 cm					
AMD 015 020 001	AMD 015 020 151	1.5	20	4F	7	14
AMD 015 020 002	AMD 015 020 152	1.5	20	4F	7	14
AMD 020 040 002	AMD 020 040 152	2.0	40	4F	7	15
AMD 020 080 002	AMD 020 080 152	2.0	80	4F	7	14
AMD 020 120 002	AMD 020 120 152	2.0	120	4F	7	14
AMD 020 150 002	AMD 020 150 152	2.0	150	4F	7	14
AMD 025 040 002	AMD 025 040 152	2.5	40	4F	7	16
AMD 025 080 002	AMD 025 080 152	2.5	80	4F	7	15
AMD 025 120 002	AMD 025 120 152	2.5	120	4F	7	14
AMD 025 150 002	AMD 025 150 152	2.5	150	4F	7	14
AMD 030 040 002	AMD 030 040 152	3.0	40	4F	7	16
AMD 030 080 002	AMD 030 080 152	3.0	80	4F	7	15
AMD 030 120 002	AMD 030 120 152	3.0	120	4F	7	14
AMD 030 150 002	AMD 030 150 152	3.0	150	4F	7	14
AMD 035 040 002	AMD 035 040 152	3.5	40	4F	7	16
AMD 035 080 002	AMD 035 080 152	3.5	80	4F	7	15
AMD 035 120 002	AMD 035 120 152	3.5	120	4F	7	14
AMD 035 150 002	AMD 035 150 152	3.5	150	4F	7	14
AMD 040 040 002	AMD 040 040 152	4.0	40	4F	7	16
AMD 040 080 002	AMD 040 080 152	4.0	80	4F	7	15
AMD 040 120 002	AMD 040 120 152	4.0	120	4F	7	14
AMD 040 150 002	AMD 040 150 152	4.0	150	4F	7	14
AMD 225 210 002	AMD 225 210 152	2.0/2.5	210	4F	7	14
AMD 253 210 002	AMD 253 210 152	2.5/3.0	210	4F	7	14
AMD 335 210 002	AMD 335 210 152	3.0/3.5	210	4F	7	14
AMD 354 210 002	AMD 354 210 152	3.5/4.0	210	4F	7	14

