

Protégé™ Rx™

Carotid Self-Expanding Stent System

The Protégé™ RX™ stent is the next generation stent designed for the anatomy of the carotid artery.

Protégé™ RX™ provides control and accurate placement for carotid interventions.

Predictable Deployment

- Proprietary EX.P.R.T.™ release technology essentially eliminates premature deployment or jumping
- No stent shortening
- Unique anatomically designed tapered stent for better fit in the carotid bifurcation
- 0.014" rapid exchange catheter with 6 F low crossing profile and flexible atraumatic tip
- Radiopaque marker on catheter clearly indicates tapered location for precise positioning

Visible Results

- Tantalum GPS™ markers enhance visibility for precise positioning and result confirmation
- Cell design produces expansion force that resists compression while providing excellent wall apposition
- Straight and tapered options for customized fit in carotid vessels

Protégé RX 6 Fr/0.014" Catheter Length 135cm

Each system includes:

One stent and delivery catheter system



Protégé™ Rx™

Carotid Self-Expanding Stent System

ORDER INFORMATION

Product Catalogue Number	Stent Dimensions		Recommended Sheath Size (F)	Recommended Guidewire (inch)	Crossing Profile (inch)
Catheter Length 135cm	Diameter (mm)	Length (mm)			
TAPERED					
SEPX-8-6-30-135	8x6	30	6	0.014	0.078
SEPX-8-6-40-135	8x6	40	6	0.014	0.078
SEPX-10-7-30-135	10x7	30	6	0.014	0.078
SEPX-10-7-40-135	10x7	40	6	0.014	0.078
STRAIGHT					
SEPX-6-20-135	6	20	6	0.014	0.078
SEPX-7-20-135	7	20	6	0.014	0.078
SEPX-8-20-135	8	20	6	0.014	0.078
SEPX-9-20-135	9	20	6	0.014	0.078
SEPX-10-20-135	10	20	6	0.014	0.078
SEPX-6-30-135	6	30	6	0.014	0.078
SEPX-7-30-135	7	30	6	0.014	0.078
SEPX-8-30-135	8	30	6	0.014	0.078
SEPX-9-30-135	9	30	6	0.014	0.078
SEPX-10-30-135	10	30	6	0.014	0.078
SEPX-6-40-135	6	40	6	0.014	0.078
SEPX-7-40-135	7	40	6	0.014	0.078
SEPX-8-40-135	8	40	6	0.014	0.078
SEPX-9-40-135	9	40	6	0.014	0.078
SEPX-10-40-135	10	40	6	0.014	0.078
SEPX-6-60-135	6	60	6	0.014	0.078
SEPX-7-60-135	7	60	6	0.014	0.078
SEPX-8-60-135	8	60	6	0.014	0.078
SEPX-9-60-135	9	60	6	0.014	0.078
SEPX-10-60-135	10	60	6	0.014	0.078

INDICATIONS: The Protégé™ RX™ is indicated for use in the iliac or subclavian arteries in the palliative treatment of malignant neoplasms in the biliary tree. It is also indicated for treatment of stenoses of the common carotid artery (CCA), internal carotid artery (ICA) and carotid bifurcation.