

## Balloon-Expandable Peripheral Stent System

- Broad offering of 6 Fr-compatible 0.035" balloon-expandable stent with radiopaque marker technology for optimized visibility.
- Low crossing profile.
- Minimal shortening for placement confidence.

VisiPro™ catheter lengths  
80 cm and 135 cm

**Each system includes:**

One stent and delivery catheter system



### COMPLIANCE CHART

Diameter (mm)	Inflation Pressure (atm)				
	8	9	10	11	12
5.0	5.00 <sup>(1)</sup>	5.09	5.16	5.22	5.28 <sup>(2)</sup>
6.0	6.00 <sup>(1)</sup>	6.11	6.22	6.31	6.39 <sup>(2)</sup>
7.0			7.00 <sup>(1)</sup>	7.09	7.17 <sup>(2)</sup>
8.0			8.00 <sup>(1)</sup>	8.15	8.26 <sup>(2)</sup>
9.0			9.00 <sup>(1)</sup>	9.15	9.28 <sup>(2)</sup>
10.0			10.00 <sup>(1)</sup>	10.11	10.21 <sup>(2)</sup>

<sup>(1)</sup> Diameter at Nominal Pressure

<sup>(2)</sup> Diameter at Rated Burst Pressure

Balloon-Expandable Peripheral Stent System

ORDER INFORMATION

Product Catalogue Number		Stent dimensions		Balloon Length (mm)	Recommended Sheath Size (F)	Recommended Guidewire (inch)	Crossing Profile (inch)
Catheter Length 80 cm	Catheter Length 135 cm	Diameter (mm)	Length (mm)				
PXP35-05-12-080							
PXP35-05-17-080	PXP35-05-17-135	5.0	17	20	6	0.035	0.079
PXP35-05-27-080	PXP35-05-27-135	5.0	27	30	6	0.035	0.079
PXP35-05-37-080	PXP35-05-37-135	5.0	37	40	6	0.035	0.079
PXP35-05-57-080	PXP35-05-57-135	5.0	57	60	6	0.035	0.079
PXP35-06-12-080		6.0	12	15	6	0.035	0.079
PXP35-06-17-080	PXP35-06-17-135	6.0	17	20	6	0.035	0.079
PXP35-06-27-080	PXP35-06-27-135	6.0	27	30	6	0.035	0.079
PXP35-06-37-080	PXP35-06-37-135	6.0	37	40	6	0.035	0.081
PXP35-06-57-080	PXP35-06-57-135	6.0	57	60	6	0.035	0.083
PXP35-07-12-080		7.0	12	15	6	0.035	0.079
PXP35-07-17-080	PXP35-07-17-135	7.0	17	20	6	0.035	0.079
PXP35-07-27-080	PXP35-07-27-135	7.0	27	30	6	0.035	0.079
PXP35-07-37-080	PXP35-07-37-135	7.0	37	40	6	0.035	0.081
PXP35-07-57-080	PXP35-07-57-135	7.0	57	60	6	0.035	0.083
PXP35-08-17-080	PXP35-08-17-135	8.0	17	20	6	0.035	0.083
PXP35-08-27-080	PXP35-08-27-135	8.0	27	30	6	0.035	0.083
PXP35-08-37-080	PXP35-08-37-135	8.0	37	40	6	0.035	0.083
PXP35-08-57-080	PXP35-08-57-135	8.0	57	60	6	0.035	0.084
PXP35-09-17-080	PXP35-09-17-135	9.0	17	20	7	0.035	0.088
PXP35-09-27-080	PXP35-09-27-135	9.0	27	30	7	0.035	0.088
PXP35-09-37-080	PXP35-09-37-135	9.0	37	40	7	0.035	0.088
PXP35-09-57-080	PXP35-09-57-135	9.0	57	60	7	0.035	0.088
PXP35-10-17-080	PXP35-10-17-135	10.0	17	20	7	0.035	0.092
PXP35-10-27-080	PXP35-10-27-135	10.0	27	30	7	0.035	0.092
PXP35-10-37-080	PXP35-10-37-135	10.0	37	40	7	0.035	0.092
PXP35-10-57-080	PXP35-10-57-135	10.0	57	60	7	0.035	0.092

Specifications Nominal  
Balloon Expandable Peripheral Stent System is indicated for use in the iliac, renal or subclavian arteries, as well as malignant biliary use.