


PROCEDURAL EFFICIENCY. OPTIMIZED.

We've taken the efficiency of our multifunctional LigaSure™ Maryland jaw device to the next level — by putting nonstick nano-coating on the jaws^{1,†}

LigaSure™ Maryland Jaw Device
with Nano-Coated Jaws



[†]Cleaning effectiveness assessed after each of two cleaning cycles.

RELIABLE VESSEL-SEALING TECHNOLOGY. NANO-COATED.

Our Maryland jaw device delivers the control and consistency of LigaSure™ vessel-sealing technology. In your hands, it gives you the power to grasp, dissect, cut, and seal.

And now with nonstick nano-coated jaws, the LigaSure™ Maryland jaw device can bring even greater efficiency to your OR.¹ Here's how.



Less sticking^{2,†}

97% less
than our legacy device

97% less
than the EnSeal™ G2 5 mm device

97% less
than the Voyant™ 5 mm Fusion

More effective cleaning

The nano-coated jaws of the LigaSure™ Maryland jaw device require fewer cleanings during a procedure compared to our legacy device. The nano-coating enables more effective jaw cleaning during the procedure than the legacy device.^{1,†}

Less eschar buildup^{2,§}

Nano-coating reduces eschar buildup by 15% compared to our legacy device.^{2,§}

Standard reprocessing methods will degrade the nonstick coating and lead to an increase in sticking, eschar buildup, and cleanings.³

Let's improve efficiency and procedural flow in your OR.

Call your local Medtronic sales representative or visit medtronic.com/covidien



1. Based on internal test report #RE00071599, LF19XX MJC marketing claims report conducted on porcine tissue. Feb. 7 to Feb. 22, 2017.

2. Based on internal test report #RE00073194, Tissue sticking comparison of the Ethicon G2™, Voyant™ 5 mm Fusion, LigaSure™ LF1737, and LigaSure™ LF1937 devices conducted on porcine tissue using the Force Triad™ energy platform. Jan. 18, 2017.

3. Based on internal report #RE00045785, LF19XX Waiver of Sticking Testing Post-reprocessing. Jan. 19, 2017.

†Cleaning effectiveness assessed after each of two cleaning cycles.

‡Instances of tissue sticking to jaws measured over 110 seals per device on the Force Triad™ energy platform.

§Eschar buildup assessed using optical imaging analysis after 60 seal and divide cycles.

© 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. ™™ Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 07/2019 – US170901(1) – [WF# 1969634]