


CELLO™
BALLOON
GUIDE
CATHETER.
EXPANDING
OPTIONS.




Medtronic
Further, Together

PRODUCT ORDERING INFORMATION

Catalog Number	Product Name	Conformable Sheath (F)	Tip Length (mm)	Balloon Length (mm)	Outer Diameter (inch)	Inner Diameter (inch)	Effective Length (cm)	Total Length (cm)
1610060	Cello 6F+	7	3	7	2.0 mm (0.079)	0.051	95	103
1610070	Cello 7F+	8	3	7	2.4 mm (0.094)	0.067	95	103
1610080	Cello 8F	8	3	7	2.7 mm (0.106)	0.075	95	103
1610090	Cello 9F	9	3	10	2.9 mm (0.114)	0.085	92	100

DEVICE DETAILS



Medtronic

Europe

Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

www.medtronic.eu

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. All claims and descriptions are for CE regulated countries. Availability of these products may vary in countries outside EU. The Cello™ balloon guide catheter is intended to temporarily block blood flow by expanding a balloon inside blood vessels during operations, such as: urgent hemostasis, hemostasis for surgery, perfusion of blood to peripheral vessel, and arterial injection for chemotherapy. *For complete instructions for use, please see IFU.

UC201706945 EE © Medtronic 2016. All rights reserved. Printed in Europe.