

# EverFlex™

## Self-Expanding Stent System

The EverFlex™ self-expanding peripheral stent system is a self-expanding nitinol stent system.

The spiral-cell interconnecting design significantly improves flexibility and vessel conformability, without sacrificing radial strength. Excellent wall apposition and compression resistance is provided by the three-wave peak design.

### **Spiral cell connection**

- Peak-to-peak connection nodes to disperse force uniformly

### **Flexible design**

- Improves fracture resistance and restores vessel patency

### **Three-wave peak design**

- Designed to resist compression and provide wall apposition

EverFlex™ catheter lengths 80cm and 120cm



## Self-Expanding Stent System

### ORDER INFORMATION

Product Catalogue Number		Stent Dimensions		Recommended		Recommended Guidewire (inch)	Crossing Profile (inch)
Catheter Length 80 cm	Catheter Length 120 cm	Diameter (mm)	Length (mm)	Lumen Size (mm)	Sheath Size (F)		
PRP35-05-020-080	PRP35-05-020-120	5	20	3.5 - 4.5	6	0.035	0.079
PRP35-05-030-080	PRP35-05-030-120	5	30	3.5 - 4.5	6	0.035	0.079
PRP35-05-040-080	PRP35-05-040-120	5	40	3.5 - 4.5	6	0.035	0.079
PRP35-05-060-080	PRP35-05-060-120	5	60	3.5 - 4.5	6	0.035	0.079
PRP35-05-080-080	PRP35-05-080-120	5	80	3.5 - 4.5	6	0.035	0.079
PRP35-05-100-080	PRP35-05-100-120	5	100	3.5 - 4.5	6	0.035	0.079
PRP35-05-120-080	PRP35-05-120-120	5	120	3.5 - 4.5	6	0.035	0.079
PRP35-05-150-080	PRP35-05-150-120	5	150	4.5 - 5.5	6	0.035	0.079
PRP35-06-020-080	PRP35-06-020-120	6	20	4.5 - 5.5	6	0.035	0.079
PRP35-06-030-080	PRP35-06-030-120	6	30	4.5 - 5.5	6	0.035	0.079
PRP35-06-040-080	PRP35-06-040-120	6	40	4.5 - 5.5	6	0.035	0.079
PRP35-06-060-080	PRP35-06-060-120	6	60	4.5 - 5.5	6	0.035	0.079
PRP35-06-080-080	PRP35-06-080-120	6	80	4.5 - 5.5	6	0.035	0.079
PRP35-06-100-080	PRP35-06-100-120	6	100	4.5 - 5.5	6	0.035	0.079
PRP35-06-120-080	PRP35-06-120-120	6	120	4.5 - 5.5	6	0.035	0.079
PRP35-06-150-080	PRP35-06-150-120	6	150	4.5 - 5.5	6	0.035	0.079
-	PRP35DR-06-200-120	6	200	4.5 - 5.5	6	0.035	0.079
PRP35-07-020-080	PRP35-07-020-120	7	20	5.5 - 6.5	6	0.035	0.079
PRP35-07-030-080	PRP35-07-030-120	7	30	5.5 - 6.5	6	0.035	0.079
PRP35-07-040-080	PRP35-07-040-120	7	40	5.5 - 6.5	6	0.035	0.079
PRP35-07-060-080	PRP35-07-060-120	7	60	5.5 - 6.5	6	0.035	0.079
PRP35-07-080-080	PRP35-07-080-120	7	80	5.5 - 6.5	6	0.035	0.079
PRP35-07-100-080	PRP35-07-100-120	7	100	5.5 - 6.5	6	0.035	0.079
PRP35-07-120-080	PRP35-07-120-120	7	120	5.5 - 6.5	6	0.035	0.079
PRP35-07-150-080	PRP35-07-150-120	7	150	5.5 - 6.5	6	0.035	0.079
-	PRP35DR-07-200-120	7	200	5.5 - 6.5	6	0.035	0.079
PRP35-08-020-080	PRP35-08-020-120	8	20	6.5 - 7.5	6	0.035	0.079
PRP35-08-030-080	PRP35-08-030-120	8	30	6.5 - 7.5	6	0.035	0.079
PRP35-08-040-080	PRP35-08-040-120	8	40	6.5 - 7.5	6	0.035	0.079
PRP35-08-060-080	PRP35-08-060-120	8	60	6.5 - 7.5	6	0.035	0.079
PRP35-08-080-080	PRP35-08-080-120	8	80	6.5 - 7.5	6	0.035	0.079
PRP35-08-100-080	PRP35-08-100-120	8	100	6.5 - 7.5	6	0.035	0.079
PRP35-08-120-080	PRP35-08-120-120	8	120	6.5 - 7.5	6	0.035	0.079
PRP35-08-150-080	PRP35-08-150-120	8	150	6.5 - 7.5	6	0.035	0.079
-	PRP35DR-08-200-120	8	200	6.5 - 7.5	6	0.035	0.079

INDICATIONS: The EverFlex™ self-expanding peripheral stent system is indicated for use in common iliac, external iliac, superficial femoral, proximal popliteal, and subclavian arteries.