


IN.PACT™ Admiral™

Paclitaxel-eluting PTA Balloon Catheter 0.035"



TECHNICAL SPECIFICATIONS

| | |
|---------------------------------|---|
| Catheter design | Over the Wire (OTW) |
| Balloon coating | FreePac – Paclitaxel and Urea (Excipient) |
| Usable shaft lengths | 40, 80 and 130 cm |
| Introducer sheath compatibility | 5, 6, 7 and 9F |
| Max. recommended guidewire | 0.035" depending on balloon size |
| Nominal pressure | 8 atm |

IN.PACT™ Admiral™

Paclitaxel-eluting PTA Balloon Catheter 0.035"

DURABLE

IN.PACT™ Admiral™ drug-eluting balloon demonstrates best-in-class clinical outcomes with durable performance through 3 years.

CONSISTENT

IN.PACT™ Admiral™ drug-eluting balloon demonstrates positive, consistent outcomes across trials, complex patients and lesion subgroups.


SAFE

IN.PACT™ Admiral™ drug-eluting balloon has an excellent safety profile, with superior results relative to PTA.

1 Balloon Diameter

2 Balloon Length

3 Usable Length



IN.PACT™ Admiral™

Paclitaxel-eluting PTA Balloon Catheter 0.035"

ORDER INFORMATION

| Product Catalogue Number | | | Balloon Diameter (mm) | Balloon Length (mm) | Recommended Introducer Sheath (F) | RBP (atm) |
|--------------------------|---------------------|----------------------|-----------------------|---------------------|-----------------------------------|-----------|
| Usable Length 40 cm | Usable Length 80 cm | Usable Length 130 cm | | | | |
| SBI 040 040 04P | SBI 040 040 08P | SBI 040 040 13P | 4 | 40 | 5 | 14 |
| SBI 040 060 04P | SBI 040 060 08P | SBI 040 060 13P | 4 | 60 | 5 | 14 |
| SBI 040 080 04P | SBI 040 080 08P | SBI 040 080 13P | 4 | 80 | 5 | 14 |
| - | SBI 040 120 08P | SBI 040 120 13P | 4 | 120 | 5 | 14 |
| - | SBI 040 150 08P | SBI 040 150 13P | 4 | 150 | 5 | 14 |
| SBI 050 040 04P | SBI 050 040 08P | SBI 050 040 13P | 5 | 40 | 6 | 14 |
| SBI 050 060 04P | SBI 050 060 08P | SBI 050 060 13P | 5 | 60 | 6 | 14 |
| SBI 050 080 04P | SBI 050 080 08P | SBI 050 080 13P | 5 | 80 | 6 | 14 |
| - | SBI 050 120 08P | SBI 050 120 13P | 5 | 120 | 6 | 14 |
| - | SBI 050 150 08P | SBI 050 150 13P | 5 | 150 | 6 | 14 |
| SBI 060 040 04P | SBI 060 040 08P | SBI 060 040 13P | 6 | 40 | 6 | 14 |
| SBI 060 060 04P | SBI 060 060 08P | SBI 060 060 13P | 6 | 60 | 6 | 14 |
| SBI 060 080 04P | SBI 060 080 08P | SBI 060 080 13P | 6 | 80 | 6 | 14 |
| - | SBI 060 120 08P | SBI 060 120 13P | 6 | 120 | 6 | 14 |
| - | SBI 060 150 08P | SBI 060 150 13P | 6 | 150 | 6 | 14 |
| SBI 070 040 04P | SBI 070 040 08P | SBI 070 040 13P | 7 | 40 | 7 | 14 |
| SBI 070 060 04P | SBI 070 060 08P | SBI 070 060 13P | 7 | 60 | 7 | 14 |
| SBI 070 080 04P | SBI 070 080 08P | SBI 070 080 13P | 7 | 80 | 7 | 14 |
| SBI 080 040 04P | SBI 080 040 08P | SBI 080 040 13P | 8 | 40 | 7 | 10 |
| SBI 080 060 04P | SBI 080 060 08P | SBI 080 060 13P | 8 | 60 | 7 | 10 |
| SBI 080 080 04P | SBI 080 080 08P | SBI 080 080 13P | 8 | 80 | 7 | 10 |
| SBI 090 040 04P | SBI 090 040 08P | SBI 090 040 13P | 9 | 40 | 7 | 10 |
| SBI 090 060 04P | SBI 090 060 08P | SBI 090 060 13P | 9 | 60 | 7 | 10 |
| SBI 090 080 04P | SBI 090 080 08P | SBI 090 080 13P | 9 | 80 | 7 | 10 |
| SBI 100 040 04P | SBI 100 040 08P | SBI 100 040 13P | 10 | 40 | 7 | 9 |
| SBI 120 040 04P | SBI 120 040 08P | SBI 120 040 13P | 12 | 40 | 9 | 9 |

IN.PACT™ Admiral™

Paclitaxel-eluting PTA Balloon Catheter 0.035"

