

IN.PACT™ Pacific™

Paclitaxel-eluting PTA Balloon Catheter 0.018"




UNPARALLELED IN PERIPHERAL

TECHNICAL SPECIFICATIONS


Catheter design	Over the Wire (OTW)
Balloon coating	FreePac – Paclitaxel and Urea (Excipient)
Usable shaft lengths	90 and 130 cm
Introducer sheath compatibility	5 F - 6 F depending on balloon size
Max. recommended guidewire	0.018"
Nominal pressure	7 atm

IN.PACT™ Pacific™

Paclitaxel-eluting PTA Balloon Catheter 0.018"



- 1** Balloon Diameter **2** Balloon Length **3** Usable Length



Paclitaxel-eluting PTA Balloon Catheter 0.018"

ORDER INFORMATION

Product Catalogue Number		Balloon Diameter (mm)	Balloon Length (mm)	Recommended Introducer Sheath (F)	RBP (atm)
Usable Length 90 cm	Usable Length 130 cm				
PCF 040 040 09P	PCF 040 040 13P	4	40	5	20
PCF 040 060 09P	PCF 040 060 13P	4	60	5	14
PCF 040 080 09P	PCF 040 080 13P	4	80	5	14
PCF 040 120 09P	PCF 040 120 13P	4	120	5	14
PCF 050 040 09P	PCF 050 040 13P	5	40	5	20
PCF 050 060 09P	PCF 050 060 13P	5	60	5	14
PCF 050 080 09P	PCF 050 080 13P	5	80	5	14
PCF 050 120 09P	PCF 050 120 13P	5	120	5	14
PCF 060 040 09P	PCF 060 040 13P	6	40	5	16
PCF 060 060 09P	PCF 060 060 13P	6	60	5	14
PCF 060 080 09P	PCF 060 080 13P	6	80	5	14
PCF 060 120 09P	PCF 060 120 13P	6	120	5	14
PCF 070 040 09P	PCF 070 040 13P	7	40	6	12
PCF 070 060 09P	PCF 070 060 13P	7	60	6	12
PCF 070 080 09P	PCF 070 080 13P	7	80	6	12
PCF 070 120 09P	PCF 070 120 13P	7	120	6	12