

Protégé™ GPS™

Self-Expanding Stent System

The Protégé™ GPS™ stent system gives control for precise stent placement avoiding jumping of the stent through the EX.P.R.T.™ retention system.

Compact delivery

- Diameters up to 14 mm
- Full line is 6 F compatible

Precision

- Proprietary EX.P.R.T.™ deployment system secures the stent to eliminate premature deployment or “jumping”
- Tantalum GPS markers enhance visibility for easier, precise positioning


Radial strength and flexibility

- Designed for radial strength without sacrificing flexibility

Each system includes:

One stent and delivery catheter system

Protégé GPS catheter lengths 80cm and 120cm



Protégé™ GPS™

Self-Expanding Stent System

ORDER INFORMATION

Product Catalogue Number		Stent Dimensions		Recommended		Recommended Guidewire (inch)	Crossing Profile (inch)
Catheter Length 80 cm	Catheter Length 120 cm	Diameter (mm)	Length (mm)	Lumen Size (mm)	Sheath Size (F)		
SERP65-09-20-80	SERP65-09-20-120	9	20	7.5 - 8.5	6	0.035	0.079
SERP65-09-30-80	SERP65-09-30-120	9	30	7.5 - 8.5	6	0.035	0.079
SERP65-09-40-80	SERP65-09-40-120	9	40	7.5 - 8.5	6	0.035	0.079
SERP65-09-60-80	SERP65-09-60-120	9	60	7.5 - 8.5	6	0.035	0.079
SERP65-09-80-80	SERP65-09-80-120	9	80	7.5 - 8.5	6	0.035	0.079
SERP65-10-20-80	SERP65-10-20-120	10	20	8.5 - 9.5	6	0.035	0.079
SERP65-10-30-80	SERP65-10-30-120	10	30	8.5 - 9.5	6	0.035	0.079
SERP65-10-40-80	SERP65-10-40-120	10	40	8.5 - 9.5	6	0.035	0.079
SERP65-10-60-80	SERP65-10-60-120	10	60	8.5 - 9.5	6	0.035	0.079
SERP65-10-80-80	SERP65-10-80-120	10	80	8.5 - 9.5	6	0.035	0.079
SERP65-12-20-80	SERP65-12-20-120	12	20	9.5 - 11.0	6	0.035	0.079
SERP65-12-30-80	SERP65-12-30-120	12	30	9.5 - 11.0	6	0.035	0.079
SERP65-12-40-80	SERP65-12-40-120	12	40	9.5 - 11.0	6	0.035	0.079
SERP65-12-60-80	SERP65-12-60-120	12	60	9.5 - 11.0	6	0.035	0.079
SERP65-12-80-80	SERP65-12-80-120	12	80	9.5 - 11.0	6	0.035	0.079
SERP65-14-20-80	SERP65-14-20-120	14	20	11.5 - 13.0	6	0.035	0.079
SERP65-14-30-80	SERP65-14-30-120	14	30	11.5 - 13.0	6	0.035	0.079
SERP65-14-40-80	SERP65-14-40-120	14	40	11.5 - 13.0	6	0.035	0.079
SERP65-14-60-80	SERP65-14-60-120	14	60	11.5 - 13.0	6	0.035	0.079
SERP65-14-80-80	SERP65-14-80-120	14	80	11.5 - 13.0	6	0.035	0.079

Specifications Nominal
INDICATIONS: The Protégé™ GPS™ stent is indicated for use in the iliac or subclavian arteries and malignant biliary use.

Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device.
Protected under one or more of the following: US Patent 6,814,746; 6,749,627; 6,623,518; 6,623,491; 6,558,415; 6,358,274; 6,132,460; D458,679. Non-US Patent pending.
Product availability and/or specifications subject to change.